

<b>IHE Integration Statement</b>		<b>Date:</b>	Dec 2011
<b>Vendor</b>	<b>Product Name</b>	<b>Version</b>	
Carestream Health	<i>Kodak DirectView</i> DR Systems (DR 9500, DR 3500, DR 7500)  Carestream DRX-1 System, DRX-Evolution, DRX-Mobile Retrofit Kit, DRX-Transportable, DRX-Ascend	5.5	
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
<b>Integration Profiles implemented</b>	<b>Actors Implemented</b>	<b>Options Implemented</b>	
Scheduled Workflow (SWF)	Acquisition Modality	None	
Patient Information Reconciliation (PIR)	Acquisition Modality	None	
Consistent Presentation of Images (CPI)	Acquisition Modality	None	
Portable Data for Imaging (PDI)	Portable Media Creator (PMC)	None	
<b>Internet address for Carestream Health IHE information:</b> <a href="http://www.carestream.com/ihe.html">http://www.carestream.com/ihe.html</a>			
<b>Links to Standards Conformance Statements for the Implementation</b>			
<b>DICOM</b>	<a href="http://www.carestream.com/dicom.html">http://www.carestream.com/dicom.html</a>		
<b>Links to general information on IHE</b>			
North America: <a href="http://www.rsna.org/IHE">http://www.rsna.org/IHE</a>	Europe: <a href="http://www.ihe-europe.org">http://www.ihe-europe.org</a>	Japan: <a href="http://www.jira-net.or.jp/ihe-j">http://www.jira-net.or.jp/ihe-j</a>	

End of Document

---

**9J6058**